

How to sample and ship a diagnostic specimen from humans

CHLAMYDIA PSITTACI diagnostics in humans

Responsible :

Name :	Prof. Dr. Daisy Vanrompay
Hospitable/Laboratory:	Universiteit Gent – Animal Science & Aquatic Ecologie
Adress :	Coupure links 653, Building B, Level -1 9000 Gent, Belgium
Tel. :	secretariat : 09/264.38.71
E-Mail :	Imbitech@ugent.be

Shipping :

Communication	<ul style="list-style-type: none">• Fill in request form (found on webpage)• Send request form by mail to the above e-mail address with confirmation of delivery date and transport conditions (sample(s) is/are frozen or not?)
Types of samples	<ul style="list-style-type: none">• pharyngeal, Broncho-alveolar lavage fluid (BAL), sputum• serum• <i>C. psittaci</i> isolates/strains for genotyping• air samples in ChlamyTrap1
Type of swab	aluminum rod, rayon-tipped
Transport medium	<ul style="list-style-type: none">• isolation; <i>C. psittaci</i> transport medium (TM)• PCR; DNA stabilization buffer (cfr. <i>C. trachomatis</i>)
Collection medium for bioaerosol monitoring	ChlamyTrap 1
Transport media available upon request	Please contact the secretariat (Imbitech@ugent.be) for transport medium or it's composition
Aerosol collection medium available upon request	Please contact the secretariat (Imbitech@ugent.be) for transport medium or it's composition
Transport for isolation	<ul style="list-style-type: none">• if delivery within 12 hours, do not freeze the samples, store at 4 °C and transport at 4 °C (ice packs).• if transport more than 12 hours, freeze (80 °C) and deliver frozen on dry ice (-70 °C) within 5 days.

Transport for PCR	4 °C following sampling, cooled transport (4°C) if possible.
Te gebruiken verpakking	<ul style="list-style-type: none"> • Conform WIV-ISP regulations • <i>C. psittaci</i> isolates/strains; biocontainer for biorisk class 3 micro-organisms or ‘triple packaging’ (IATA 602 instructions)
Maximum duration of transport	<ul style="list-style-type: none"> • for isolation 12 hours (4 °C) to max 5 days (-70 °C) • for PCR < 1 week

Tests :

Type of test	routine test	upon request
Confirmation of diagnosis	X	
Diagnosis (culture)	X	
Diagnosis (PCR)	X	
Molecular identification		X
Serologic diagnosis	X	
Billing of these tests to the applying laboratory	X	X

Results :

Number of days until result of a routine test	<ul style="list-style-type: none"> • Culture: 1 to 2 weeks • PCR: 2 to 4 days
Number of days until result of an analysis upon request	1 to 4 weeks